

RESMED

**Lasting
Results**

QUALITY SYSTEMS

COMPLETE SOLUTIONS

RELIABLE PRODUCTS

**ResMed Corp
Product List**

**Effective: January 2012
Contents**

- Swift FX Nasal Pillows System and Accessories.....2
- Swift FX for Her and Swift FX Bella Nasal Pillows Systems.....3
- Swift LT and Swift LT for Her Nasal Pillows Systems.....4,5
- Mirage Swift II Nasal Pillows System and Accessories.....6
- Mirage FX and Mirage FX for Her Nasal Mask Systems.....7
- Mirage SoftGel Nasal Mask System and Accessories.....8,9
- Mirage Activa LT Nasal Mask System and Accessories.....10,11
- Mirage SoftGel/Activa LT ConvertAble Pack.....12
- Mirage Activa Nasal Mask System and Accessories.....13
- Mirage Micro and Mirage Micro for Kids Nasal Mask Systems.....14
- Ultra Mirage II Nasal Mask System and Accessories.....15
- Pixi Pediatric Nasal Mask System.....16
- Mirage Vista and Mirage Kidsta Nasal Mask Systems.....17
- Quattro FX and Quattro FX for Her Full Face Mask Systems....18,19
- Mirage Quattro Full Face Mask System and Accessories.....20,21
- Ultra Mirage Full Face Mask System and Accessories.....22,23
- Mirage Liberty Full Face Mask System and Accessories.....24,25
- Mask Spares and Accessories.....26
- S9 APAP and CPAP Flow Generators.....27,28
- S9 VPAP Bilevel Flow Generators.....29
- S9 Humidifiers and Accessories.....30
- Ventilation Products.....31
- C-Series Flow Generator.....32
- Software, Modules and Accessories.....33,34
- ApneaLink/ApneaLink Plus and Accessories.....35,36
- Flow Generator Accessories and Filters.....37,38,39
- Heated Humidification Systems and Accessories.....40



Item #	Description	HCPCS
61500	Swift™ FX nasal pillows system Includes small, medium, and large size dual-wall pillows, and headgear	A7034, +A7035
PILLOWS		
61520	Extra small dual-wall pillow	A7033
61521	Small dual-wall pillow	
61522	Medium dual-wall pillow	
61523	Large dual-wall pillow	
FRAME		
61510	Swift FX frame system with extra small dual-wall pillows - no headgear	A7034
61511	Swift FX frame system with small dual-wall pillows - no headgear	
61512	Swift FX frame system with medium dual-wall pillows - no headgear	
61513	Swift FX frame system with large dual-wall pillows - no headgear	
61528	Swift FX frame assembly - no dual-wall pillows, no headgear	N/A
HEADGEAR		
61529	Headgear assembly	A7035
ACCESSORIES		
61528	Short tube assembly - refer to Swift FX frame assembly (Item #61528)	N/A
61530	Soft wraps - 2 pack	N/A
61121	Distal pressure port attachment	N/A
60549	Pressure port plug - 1 pack	N/A
60550	Pressure port plug - 10 pack	N/A

**Swift™ FX for Her
Swift™ FX Bella
Nasal Pillows Systems**



Item #	Description	HCPCS
61540	Swift™ FX for Her nasal pillows system Includes extra small, small, and medium size dual-wall pillows, and headgear	A7034, +A7035
61560	Swift™ FX Bella nasal pillows system Includes extra small, small, and medium size dual-wall pillows, Bella Loops accessory and Swift FX for Her headgear with soft wraps	
PILLOWS		
61520	Extra small dual-wall pillow	A7033
61521	Small dual-wall pillow	
61522	Medium dual-wall pillow	
61523	Large dual-wall pillow	
FRAME		
61510	Swift FX frame system with extra small dual-wall pillows - no headgear	A7034
61511	Swift FX frame system with small dual-wall pillows - no headgear	
61512	Swift FX frame system with medium dual-wall pillows - no headgear	
61513	Swift FX frame system with large dual-wall pillows - no headgear	
61528	Swift FX frame assembly - no dual-wall pillows, no headgear	N/A
HEADGEAR		
61543	Swift FX for Her headgear assembly (headgear, pink and light blue back straps, pink and gray soft wraps)	A7035
61580	Swift FX for Her headgear assembly with Bella Loops accessory (headgear, pink and light blue back straps, pink and gray soft wraps, Bella Loops)	
ACCESSORIES		
61528	Short tube assembly - refer to Swift FX frame assembly (Item #61528)	N/A
61530	Soft wraps - 2 pack (gray)	N/A
61544	Soft wraps - 2 pack (pink)	N/A
61121	Distal pressure port attachment	N/A
60549	Pressure port plug - 1 pack	N/A
60550	Pressure port plug - 10 pack	N/A
61581	Swift FX Bella Loops	N/A



Item #	Description	HCPCS
60560	Swift™ LT nasal pillows system Includes small, medium, and large size dual-wall pillows, and headgear	A7034, +A7035
PILLOWS		
60574	Extra small dual-wall pillow	A7033
60571	Small dual-wall pillow	
60572	Medium dual-wall pillow	
60573	Large dual-wall pillow	
FRAME		
60585	Swift LT frame system with extra small dual-wall pillows - no headgear	A7034
60582	Swift LT frame system with small dual-wall pillows - no headgear	
60583	Swift LT frame system with medium dual-wall pillows - no headgear	
60584	Swift LT frame system with large dual-wall pillows - no headgear	
60581	Swift LT frame assembly - no dual-wall pillows, no headgear	N/A
60575	Frame only	N/A
HEADGEAR		
60578	Headgear assembly	A7035
ACCESSORIES		
61121	Distal pressure port attachment	N/A
60558	Headgear buckle (back) - 1 pack	N/A
60559	Headgear buckle (back) - 10 pack	N/A
60579	Headgear buckle (top) - 1 pack	N/A
60580	Headgear buckle (top) - 10 pack	N/A
60577	Short tube assembly, includes elbow, short tube and swivel - 1 pack	N/A
60742	Short tube assembly, includes elbow, short tube and swivel - 10 pack	N/A
60586	Tube retainer (plastic)	N/A
60557	Tube retainer (velcro)	N/A
60549	Pressure port plug - 1 pack	N/A
60550	Pressure port plug - 10 pack	N/A
60599	Soft wraps - 2 pack	N/A



Item #	Description	HCPCS
60588	Swift™ LT for Her nasal pillows system Includes extra small, small and medium size dual-wall pillows, and headgear	A7034, +A7035
PILLOWS		
60574	Extra small dual-wall pillow	A7033
60571	Small dual-wall pillow	
60572	Medium dual-wall pillow	
60573	Large dual-wall pillow	
FRAME		
60585	Swift LT frame system with extra small dual-wall pillows - no headgear	A7034
60582	Swift LT frame system with small dual-wall pillows - no headgear	
60583	Swift LT frame system with medium dual-wall pillows - no headgear	
60584	Swift LT frame system with large dual-wall pillows - no headgear	
60581	Swift LT frame assembly - no dual-wall pillows, no headgear	N/A
60575	Frame only	N/A
HEADGEAR		
60595	Headgear assembly - light blue	A7035
ACCESSORIES		
61121	Distal pressure port attachment	N/A
60740	Headgear buckle (back) - 1 pack - light blue	N/A
60741	Headgear buckle (back) - 10 pack - light blue	N/A
60597	Headgear buckle (top) - 1 pack - light blue	N/A
60596	Headgear buckle (top) - 10 pack - light blue	N/A
60577	Short tube assembly, includes elbow, short tube and swivel - 1 pack	N/A
60742	Short tube assembly, includes elbow, short tube and swivel - 10 pack	N/A
60598	Tube retainer (plastic) - light blue	N/A
60549	Pressure port plug - 1 pack	N/A
60550	Pressure port plug - 10 pack	N/A
60599	Soft wraps - 2 pack	N/A



Item #	Description	HCPCS
60512	Mirage Swift™ II nasal pillows system Includes small, medium, and large size pillow sleeves, and headgear	A7034, +A7035
PILLOW SLEEVES		
60541	Small pillow sleeve	A7033
60542	Medium pillow sleeve	
60543	Large pillow sleeve	
FRAME		
60537	Mirage Swift II frame system with small pillows - no headgear	A7034
60538	Mirage Swift II frame system with medium pillows - no headgear	
60539	Mirage Swift II frame system with large pillows - no headgear	
60544	Mirage Swift II frame assembly - no pillows, no headgear	N/A
60536	Frame with clip and cap only	N/A
HEADGEAR		
60526	Headgear	A7035
ACCESSORIES		
61121	Distal pressure port attachment	N/A
60523	Headgear buckle	N/A
60540	Frame clip	N/A
60524	Seal ring - 10 pack	N/A
60525	Seal cap	N/A
60528	Short tube assembly, includes elbow, short tube and swivel - 1 pack	N/A
60529	Short tube assembly, includes elbow, short tube and swivel - 10 pack	N/A
60530	Tube retainer - 3 pack	N/A
60549	Pressure port plug - 1 pack	N/A
60550	Pressure port plug - 10 pack	N/A



Item #	Description	HCPCS
62103	Mirage™ FX nasal mask complete system Includes cushion and headgear	
62118	Mirage FX nasal mask complete system - wide Includes cushion and headgear	A7034, +A7035
62128	Mirage™ FX for Her nasal mask complete system Includes cushion and headgear	
CUSHIONS		
62111	Mirage FX and Mirage FX for Her cushion - standard	A7032
62125	Mirage FX cushion - wide	
FRAME SYSTEM		
62112	Mirage FX and Mirage FX for Her frame system with cushion - standard	A7034
62126	Mirage FX frame system with cushion - wide	
HEADGEAR		
62110	Mirage FX headgear	A7035
62129	Mirage FX for Her headgear	



Item #	Description	HCPCS
61600	Mirage™ SoftGel nasal mask complete system - small Includes cushion and headgear	
61601	Mirage SoftGel nasal mask complete system - medium Includes cushion and headgear	A7034,
61602	Mirage SoftGel nasal mask complete system - large Includes cushion and headgear	+A7035
61603	Mirage SoftGel nasal mask complete system - large wide Includes cushion and headgear	
CUSHIONS		
61631	Small cushion and clip	
61632	Medium cushion and clip	A7032
61633	Large cushion and clip	
61634	Large wide cushion and clip	
FRAME		
61626	Mirage SoftGel nasal mask frame system with small cushion - no headgear *	
61627	Mirage SoftGel nasal mask frame system with medium cushion - no headgear *	A7034
61628	Mirage SoftGel nasal mask frame system with large cushion - no headgear *	
61629	Mirage SoftGel nasal mask frame system with large wide cushion - no headgear * * includes headgear clips	
60196	Mirage Activa LT/SoftGel nasal mask complete frame assembly - small, no cushion, no headgear, includes headgear clips	N/A
60175	Mirage Activa LT/SoftGel nasal mask complete frame assembly - standard, no cushion, no headgear, includes headgear clips	
60180	Mask frame only - standard	N/A

Item #	HEADGEAR	HCPCS
16117	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, and Quattro - medium	A7035
16118	Universal headgear for Ultra Mirage II, Micro, Micro for Kids, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - small	
16119	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro- large	A7035
ACCESSORIES		
16399	Elbow assembly - includes elbow and swivel	N/A
60176	Cushion clip - standard	N/A
60197	Cushion clip - small	N/A
60123	Forehead pad - 1 pack	N/A
60124	Forehead pad - 10 pack	N/A
16393	Forehead support	N/A
16569	Headgear clip - 2 pack	N/A
16734	Headgear clip - 10 pack	N/A
16570	Port cap - 2 pack	N/A
16571	Port cap - 10 pack	N/A
16565	Quick release swivel	N/A
61289	Forehead Dial	N/A
16387	Elbow	N/A



Item #	Description	HCPCS
60182	Mirage Activa™ LT nasal mask complete system - small Includes cushion and headgear	
60148	Mirage Activa LT nasal mask complete system - medium Includes cushion and headgear	A7034,
60149	Mirage Activa LT nasal mask complete system - large Includes cushion and headgear	+A7035
60150	Mirage Activa LT nasal mask complete system - large wide Includes cushion and headgear	
CUSHIONS		
60198	Small cushion and clip	
60177	Medium cushion and clip	A7032
60178	Large cushion and clip	
60179	Large wide cushion and clip	
FRAME		
60195	Mirage Activa LT nasal mask frame system with small cushion - no headgear *	
60172	Mirage Activa LT nasal mask frame system with medium cushion - no headgear *	
60173	Mirage Activa LT nasal mask frame system with large cushion - no headgear *	A7034
60174	Mirage Activa LT nasal mask frame system with large wide cushion - no headgear *	
	* includes headgear clips	
60196	Mirage Activa LT/SoftGel nasal mask complete frame assembly - small, no cushion, no headgear, includes headgear clips	N/A
60175	Mirage Activa LT/SoftGel nasal mask complete frame assembly - standard, no cushion, no headgear, includes headgear clips	
60193	Mask frame only - small	
60180	Mask frame only - standard	N/A

**Mirage Activa LT
Nasal Mask Systems (cont.)**

Item #	HEADGEAR	HCPCS
16117	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, and Quattro - medium	A7035
16118	Universal headgear for Ultra Mirage II, Micro, Micro for Kids, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - small	A7035
16119	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro- large	
ACCESSORIES		
16399	Elbow assembly - includes elbow and swivel	N/A
60176	Cushion clip - standard	N/A
60197	Cushion clip - small	N/A
60123	Forehead pad - 1 pack	N/A
60124	Forehead pad - 10 pack	N/A
16393	Forehead support	N/A
16569	Headgear clip - 2 pack	N/A
16734	Headgear clip - 10 pack	N/A
16570	Port cap - 2 pack	N/A
16571	Port cap - 10 pack	N/A
16565	Quick release swivel	N/A
61289	Forehead Dial	N/A
16387	Elbow	N/A

Item #	Description	HCPCS
61604	Mirage SoftGel/Activa LT ConvertAble™ Pack - small Includes one SoftGel mask, one Activa LT cushion and headgear	
61609	Mirage SoftGel/Activa LT ConvertAble Pack - medium Includes one SoftGel mask, one Activa LT cushion and headgear	A7032, +A7034, +A7035
61615	Mirage SoftGel/Activa LT ConvertAble Pack - large Includes one SoftGel mask, one Activa LT cushion and headgear	
61620	Mirage SoftGel/Activa LT ConvertAble Pack - large wide Includes one SoftGel mask, one Activa LT cushion and headgear	





Item #	Description	HCPCS
60100	Mirage Activa™ nasal mask complete system - standard Includes cushion and headgear	
60101	Mirage Activa nasal mask complete system - large Includes cushion and headgear	A7034, +A7035
60102	Mirage Activa nasal mask complete system - shallow Includes cushion and headgear	
CUSHIONS		
60117	Standard cushion	
60118	Large cushion	A7032
60119	Shallow cushion	
FRAME		
60111	Mirage Activa nasal mask frame system with standard cushion - no headgear *	A7034
60112	Mirage Activa nasal mask frame system with large cushion - no headgear *	
60113	Mirage Activa nasal mask frame system with shallow cushion - no headgear * * includes headgear clips	
60144	Mirage Activa nasal mask complete frame assembly - no cushion, no headgear, includes headgear clips	N/A
60121	Mask frame - one size fits all, mask frame only	N/A
HEADGEAR		
16117	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, and Quattro - medium	A7035
16118	Universal headgear for Ultra Mirage II, Micro, Micro for Kids, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - small	A7035
16119	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - large	
ACCESSORIES		
60126	Carry bag	N/A
60120	Cushion clip	N/A
60125	Elbow assembly - includes exhaust vent	N/A
60123	Forehead pad - 1 pack	N/A
60124	Forehead pad - 10 pack	N/A
60122	Forehead support	N/A
60115	Headgear clips - 2 pack	N/A
60116	Headgear clips - 10 pack	N/A
60926	Inlet tube with swivel assembly - 1 pack	N/A
60927	Inlet tube with swivel assembly - 10 pack	N/A
16570	Port cap - 2 pack	N/A
16571	Port cap - 10 pack	N/A
60014	Swivel assembly - 10 pack	N/A
61905	Mask fitting template	N/A



Item #	Description	HCPCS
61013	Mirage Micro™ for Kids nasal mask complete system - small Includes small cushion and small headgear	
16333	Mirage Micro™ nasal mask complete system - small Includes small cushion and medium headgear	A7034, A7035
16334	Mirage Micro nasal mask complete system - medium and large Includes two cushions (medium and large) and medium headgear	
16335	Mirage Micro nasal mask complete system - large wide and x-large Includes two cushions (large wide and x-large) and medium headgear	
CUSHIONS		
16388	Small cushion	
16389	Medium cushion	
16390	Large cushion	A7032
16391	Large wide cushion	
16392	X-large cushion	
FRAME		
16378	Mirage Micro and Mirage Micro for Kids nasal mask frame system with small cushion - no headgear	
16379	Mirage Micro nasal mask frame system with medium cushion - no headgear	A7034
16380	Mirage Micro nasal mask frame system with large cushion - no headgear	
16381	Mirage Micro nasal mask frame system with large wide cushion - no headgear	
16382	Mirage Micro nasal mask frame system with x-large cushion - no headgear	
16383	Mirage Micro and Mirage Micro for Kids nasal mask complete frame assembly, small - no cushion, no headgear	N/A
16384	Mirage Micro nasal mask complete frame assembly, standard - no cushion, no headgear	
16394	Mask frame only, standard	N/A
16395	Mask frame only, small	
HEADGEAR		
16117	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, and Quattro - medium	A7035
16118	Universal headgear for Ultra Mirage II, Micro, Micro for Kids, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - small	A7035
16119	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - large	
ACCESSORIES		
16399	Elbow assembly - includes elbow and swivel	N/A
60123	Forehead pad - 1 pack	N/A
60124	Forehead pad - 10 pack	N/A
16393	Forehead support	N/A
16569	Headgear clip - 2 pack	N/A
16734	Headgear clip - 10 pack	N/A
16570	Port cap - 2 pack	N/A
16571	Port cap - 10 pack	N/A
16565	Quick release swivel	N/A
61289	Forehead Dial	N/A
16387	Elbow	N/A



Item #	Description	HCPCS
16548	Ultra Mirage™ II nasal mask complete system – standard Includes cushion and headgear	
16549	Ultra Mirage II nasal mask complete system – large Includes cushion and headgear	
16550	Ultra Mirage II nasal mask complete system – shallow Includes cushion and headgear	A7034, +A7035
16577	Ultra Mirage II nasal mask complete system – shallow wide Includes cushion and headgear	
CUSHIONS		
16556	Standard cushion	
16557	Large cushion	
16558	Shallow cushion	A7032
16735	Shallow wide cushion	
FRAME		
16726	Ultra Mirage II nasal mask frame system with standard cushion - no headgear	
16727	Ultra Mirage II nasal mask frame system with shallow cushion - no headgear	
16728	Ultra Mirage II nasal mask frame system with large cushion - no headgear	A7034
16736	Ultra Mirage II nasal mask frame system with shallow wide cushion - no headgear	
16594	Ultra Mirage II nasal mask complete frame assembly - no cushion, no headgear	N/A
HEADGEAR		
16117	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, and Quattro - medium	A7035
16118	Universal headgear for Ultra Mirage II, Micro, Micro for Kids, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - small	
16119	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - large	A7035
ACCESSORIES		
16563	Cushion clip	N/A
16566	Elbow assembly - includes elbow, quick release swivel and vent cover	N/A
16572	Elbow clip - 1 pack	N/A
16573	Elbow clip - 10 pack	N/A
60123	Forehead pad - 1 pack	N/A
60124	Forehead pad - 10 pack	N/A
16590	Forehead support with pads	N/A
16569	Headgear release clip - 2 pack	N/A
16734	Headgear release clip - 10 pack	N/A
16570	Port cap - 2 pack	N/A
16571	Port cap - 10 pack	N/A
16565	Quick release swivel	N/A
16575	Vent cover	N/A
61905	Mask fitting template	N/A



Item #	Description	HCPCS
61030	Pixi™ Pediatric nasal mask complete system Includes cushion and headgear	A7034, +A7035
CUSHIONS		
61033	Pixi Pediatric cushion	A7032
FRAME SYSTEM		
61035	Pixi Pediatric frame system - no headgear	A7034
HEADGEAR		
61034	Pixi Pediatric headgear	A7035



The Mirage Kidsta™ nasal mask is intended for multi-patient use for children aged seven and older or weighing more than 40 lbs, prescribed continuous positive airway pressure (CPAP) or bilevel therapy in hospital, clinic and/or home environments. The Mirage Kidsta nasal mask is intended to be used with ResMed's VPAP™ III ST-A.

Item #	Description	HCPCS
60000	Mirage Vista™ nasal mask complete system - standard Includes cushion and headgear	
60001	Mirage Vista nasal mask complete system - deep Includes cushion and headgear	A7034, +A7035
61010	Mirage Kidsta™ nasal mask complete system - small Includes cushion and headgear	
CUSHIONS		
60921	Mirage Vista Standard cushion	
60922	Mirage Vista Deep cushion	A7032
60013	Mirage Kidsta Small Cushion	
FRAME		
60929	Mirage Vista nasal mask frame system with standard cushion - no headgear	
60930	Mirage Vista nasal mask frame system with deep cushion - no headgear	A7034
60928	Mirage Vista nasal mask complete frame assembly- no cushion, no headgear	N/A
60924	Mirage Vista and Mirage Kidsta Mask frame - one size fits all, mask frame only	N/A
HEADGEAR		
60918	Mirage Vista Headgear assembly	
61020	Mirage Kidsta Headgear assembly	A7035
ACCESSORIES		
60919	Mirage Vista Headgear clips - 2 pack	N/A
60920	Mirage Vista Headgear clips - 10 pack	N/A
61025	Mirage Kidsta Headgear clips - 2 pack	N/A
61026	Mirage Kidsta Headgear clips - 10 pack	N/A
60926	Inlet tube with swivel assembly - 1 pack	N/A
60927	Inlet tube with swivel assembly - 10 pack	N/A
16570	Port caps - 2 pack	N/A
16571	Port caps - 10 pack	N/A
60014	Swivel assembly - 10 pack	N/A
60940	2D fitting template	N/A



Item #	Description	HCPCS
61700	Quattro™ FX full face mask complete system - small Includes cushion and medium/standard headgear	
61701	Quattro™ FX full face mask complete system - medium Includes cushion and medium/standard headgear	A7030, +A7035
61702	Quattro™ FX full face mask complete system - large Includes cushion and medium/standard headgear	
CUSHIONS		
61721	Small cushion	
61722	Medium cushion	
61723	Large cushion	A7031
FRAME SYSTEM		
61724	Quattro FX full face mask frame system with small cushion - no headgear	
61725	Quattro FX full face mask frame system with medium cushion - no headgear	
61726	Quattro FX full face mask frame system with large cushion - no headgear	A7030
HEADGEAR		
61733	Quattro FX full face mask headgear - medium/standard	
61734	Quattro FX full face mask headgear - small	
61738	Quattro FX full face mask headgear - large	A7035



Item #	Description	HCPCS
62501	Quattro™ FX for Her full face mask complete system - small Includes cushion and medium/standard headgear	A7030,
62502	Quattro™ FX for Her full face mask complete system - medium Includes cushion and medium/standard headgear	+A7035
CUSHIONS		
61721	Small cushion	A7031
61722	Medium cushion	
61723	Large cushion	
FRAME SYSTEM		
62513	Quattro FX for Her full face mask frame system with small cushion - no headgear	A7030
62514	Quattro FX for Her full face mask frame system with medium cushion - no headgear	
HEADGEAR		
62506	Quattro FX for Her full face mask headgear (light pink) - small	A7035
62507	Quattro FX for Her full face mask headgear (light pink) - standard	



Item #	Description	HCPCS
61200	Mirage Quattro™ full face mask complete system - x-small Includes cushion and medium headgear	
61201	Mirage Quattro full face mask complete system - small Includes cushion and medium headgear	A7030,
61202	Mirage Quattro full face mask complete system - medium Includes cushion and medium headgear	+A7035
61203	Mirage Quattro full face mask complete system - large Includes cushion and medium headgear	
CUSHIONS		
61290	X-small cushion with clip	
61291	Small cushion with clip	A7031
61292	Medium cushion with clip	
61293	Large cushion with clip	
FRAME		
61260	Mirage Quattro full face mask frame system with x-small cushion - no headgear	
61261	Mirage Quattro full face mask frame system with small cushion - no headgear	A7030
61262	Mirage Quattro full face mask frame system with medium cushion - no headgear	
61263	Mirage Quattro full face mask frame system with large cushion - no headgear	
61264	Mirage Quattro full face mask complete frame assembly, x-small - no cushion, no headgear	
61265	Mirage Quattro full face mask complete frame assembly, small - no cushion, no headgear	
61266	Mirage Quattro full face mask complete frame assembly, medium - no cushion, no headgear	N/A
61267	Mirage Quattro full face mask complete frame assembly, large - no cushion, no headgear	
61270	Mask frame only, x-small	
61271	Mask frame only, small	N/A
61272	Mask frame only, medium	
61273	Mask frame only, large	

Mirage Quattro Full Face Mask Systems (cont.)

Item #	HEADGEAR	HCPCS
16117	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, and Quattro - medium	A7035
16118	Universal headgear for Ultra Mirage II, Micro, Micro for Kids, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - small	A7035
16119	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - large	
ACCESSORIES		
61286	Valve and clip assembly - 1 pack	N/A
61287	Valve and clip assembly - 10 pack	N/A
61274	Cushion clip, x-small	N/A
61275	Cushion clip, small	N/A
61276	Cushion clip, medium	N/A
61277	Cushion clip, large	N/A
61294	Elbow assembly - includes elbow, valve, clip and swivel	N/A
61282	Elbow	N/A
60123	Forehead pad - 1 pack	N/A
60124	Forehead pad - 10 pack	N/A
61288	Forehead support	N/A
61289	Forehead dial	N/A
16565	Quick release swivel	N/A
60115	Headgear clips - 2 pack	N/A
60116	Headgear clips - 10 pack	N/A
61295	Ports cap - 2 pack	N/A
61296	Ports cap - 10 pack	N/A
61858	Fitting template	N/A



Item #	Description	HCPCS
60600	Ultra Mirage™ full face mask complete system - small standard Includes cushion and headgear	
60601	Ultra Mirage full face mask complete system - small shallow Includes cushion and headgear	
60602	Ultra Mirage full face mask complete system - medium standard Includes cushion and headgear	
60603	Ultra Mirage full face mask complete system - medium shallow Includes cushion and headgear	A7030, +A7035
60604	Ultra Mirage full face mask complete system - large standard Includes cushion and headgear	
60605	Ultra Mirage full face mask complete system - large shallow Includes cushion and headgear	
CUSHIONS		
16604	Small standard cushion	
16671	Small shallow cushion	
16605	Medium standard cushion	
16672	Medium shallow cushion	A7031
16606	Large standard cushion	
16673	Large shallow cushion	
FRAME		
60620	Ultra Mirage full face mask frame system with small standard cushion - no headgear	
60621	Ultra Mirage full face mask frame system with small shallow cushion - no headgear	
60622	Ultra Mirage full face mask frame system with medium standard cushion - no headgear	A7030
60623	Ultra Mirage full face mask frame system with medium shallow cushion - no headgear	
60624	Ultra Mirage full face mask frame system with large standard cushion - no headgear	
60625	Ultra Mirage full face mask frame system with large shallow cushion - no headgear	
60666	Ultra Mirage full face mask complete frame assembly, small - no cushion, no headgear	
60667	Ultra Mirage full face mask complete frame assembly, medium - no cushion, no headgear	N/A
60668	Ultra Mirage full face mask complete frame assembly, large - no cushion, no headgear	
60670	Mask frame only, small	
60671	Mask frame only, medium	N/A
60672	Mask frame only, large	

Item #	HEADGEAR	HCPCS
60674	Headgear - includes headgear clips, one size fits all	A7035
16118	Universal headgear for Ultra Mirage II, Micro, Micro for Kids, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - small	A7035
16119	Universal headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel, Ultra Mirage Full Face Mask and Quattro - large	
ACCESSORIES		
16613	Anti-asphyxia valve - includes elbow retainer, valve frame, valve membrane, valve elbow, swivel clip and swivel	N/A
16617	Anti-asphyxia valve membrane - 1 pack	N/A
16618	Anti-asphyxia valve membrane - 10 pack	N/A
16674	Cushion clip, small	N/A
16675	Cushion clip, medium	N/A
16676	Cushion clip, large	N/A
16637	Elbow and swivel clip kit - includes five elbow clips and five quick release clips	N/A
60683	Elbow clip	N/A
60123	Forehead pad - 1 pack	N/A
60124	Forehead pad - 10 pack	N/A
60673	Forehead support - 1 pack	N/A
60115	Headgear clips - 2 pack	N/A
60116	Headgear clips - 10 pack	N/A
60681	Luer locks - 2 pack	N/A
60682	Luer locks - 10 pack	N/A
16689	Swivel only	N/A
60680	Swivel clip	N/A
16624	Valve elbow	N/A
16981	3D fitting template	N/A





Item #	Description	HCPCS
61300	Mirage Liberty™ full face mask complete system - small Includes cushion, pillows and headgear	A7027, +A7035
61301	Mirage Liberty full face mask complete system - large Includes cushion, pillows and headgear	
CUSHIONS		
61330	Small mouth cushion	A7028
61331	Large mouth cushion	
PILLOW SLEEVES		
61333	Dual-wall nasal pillows - small	A7029
61334	Dual-wall nasal pillows - medium	
61335	Dual-wall nasal pillows - large	
FRAME		
61318	Mirage Liberty full face mask frame system with small cushion and small pillows - no headgear	A7027
61319	Mirage Liberty full face mask frame system with small cushion and medium pillows - no headgear	
61320	Mirage Liberty full face mask frame system with small cushion and large pillows - no headgear	
61321	Mirage Liberty full face mask frame system with large cushion and small pillows - no headgear	
61322	Mirage Liberty full face mask frame system with large cushion and medium pillows - no headgear	
61323	Mirage Liberty full face mask frame system with large cushion and large pillows - no headgear	
61325	Mirage Liberty full face mask complete frame assembly, small - no cushion, no pillows, no headgear	N/A
61326	Mirage Liberty full face mask complete frame assembly, large - no cushion, no pillows, no headgear	
61337	Mask frame only, small	N/A
61338	Mask frame only, large	

Mirage Liberty Full Face Mask Systems (cont.)

Item #	HEADGEAR	HCPCS
61348	Headgear - small	A7035
61349	Headgear - standard	
ACCESSORIES		
61342	Valve and clip assembly - 1 pack	N/A
61343	Valve and clip assembly - 10 pack	N/A
61340	Elbow assembly - includes elbow, valve, clip and swivel	N/A
61352	Upper clips (left & right)	N/A
61353	Lower clips - 2 pack	N/A
16570	Ports cap - 2 pack	N/A
16571	Ports cap - 10 pack	N/A
60926	Inlet tube with swivel assembly - 1 pack	N/A
60927	Inlet tube with swivel assembly - 10 pack	N/A
61900	Fitting template	N/A
60014	Swivel assembly - 10 pack	N/A

Item #	Description	HCPCS
16015	Chin Restraint	A7036
61914 61915	<p>Gecko™ Nasal Pad: Small - 10 pack (case) Gecko™ Nasal Pad: Large - 10 pack (case)</p> <p>Note: The Gecko Nasal Pad is reimbursed on the individual units: PN 61911 (small) and PN 61910 (large). The Gecko Nasal Pad can only be ordered and shipped by 10 pack (case). The individual item numbers cannot be ordered.</p>	A7032
16120	Camouflage Headgear - Medium - Universal Headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel and Quattro masks	A7035
16121	Light Blue Headgear - Medium - Universal Headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel and Quattro masks	
16122	Pink Headgear - Medium - Universal Headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel and Quattro masks	
16123	Leopard Headgear - Medium - Universal Headgear for Ultra Mirage II, Micro, Activa, Activa LT, SoftGel and Quattro masks	



The S9™ system sets the new industry standard for comfortable sleep apnea treatment. Equipped with sophisticated new features that showcase ResMed's innovation and dedication to patient compliance, the stylish S9 system delivers everything you and your patients need - and more.

Item #	Description	HCPCS
36005	<p>S9 AutoSet™ Features include the new Enhanced AutoSet algorithm with central sleep apnea detection. Superior APAP comfort with whisper quiet operation and dramatically reduced conducted noise levels. Also includes EPR with Easy-Breathe technology; available in AutoSet mode as well as CPAP mode. Includes usage-based compliance and efficacy downloads for accurate patient monitoring and a standard data card for easy two-way data transfer. Optional accessories support wired and wireless compliance data transfer and integrated pulse oximetry monitoring.</p>	E0601
36015	<p>S9 AutoSet & H5i™ S9 AutoSet packaged with the H5i humidifier (Item #36900).</p>	E0601, +E0562
36075	<p>S9 AutoSet & H5i for Her (light pink) S9 AutoSet for Her packaged with the H5i for Her humidifier.</p>	E0601, +E0562
36025	<p>S9 AutoSet & Climate Control Kit S9 AutoSet packaged with the H5i humidifier and the ClimateLine™ Tubing (Item #36995)</p>	E0601, +E0562, +A4604
36003	<p>S9 Elite™ Features include new central sleep apnea detection, whisper quiet operation and dramatically reduced conducted noise levels. Also includes EPR with Easy-Breathe technology; Includes usage-based compliance and efficacy downloads for accurate patient monitoring and a standard data card for easy two-way data transfer. Optional accessories support wired and wireless compliance data transfer and integrated pulse oximetry monitoring.</p>	E0601
36013	<p>S9 Elite & H5i S9 Elite packaged with the H5i humidifier (Item #36900).</p>	E0601, +E0562
36023	<p>S9 Elite & Climate Control Kit S9 Elite packaged with the H5i humidifier and the ClimateLine Tubing (Item #36995)</p>	E0601, +E0562, +A4604



The S9™ system sets the new industry standard for comfortable sleep apnea treatment. Equipped with sophisticated new features that showcase ResMed's innovation and dedication to patient compliance, the stylish S9 system delivers everything you and your patients need - and more.

Item #	Description	HCPCS
36002	<p>S9 Escape Auto™ Features include superior APAP mode based upon the gold-standard AutoSet Algorithm providing comfort with whisper quiet operation and dramatically reduced conducted noise levels. Also includes EPR with Easy-Breathe technology; available in AutoSet mode as well as CPAP mode. Includes usage-based compliance and AHI data for accurate patient monitoring, along with a standard data card for easy two-way data transfer. Optional accessories support wired and wireless compliance data transfer.</p>	E0601
36012	<p>S9 Escape Auto & H5i S9 Escape Auto packaged with the H5i humidifier (Item #36900).</p>	E0601, +E0562
36022	<p>S9 Escape Auto & Climate Control Kit S9 Escape Auto packaged with the H5i humidifier and the ClimateLine Tubing (Item #36995)</p>	E0601, +E0562, +A4604
36001	<p>S9 Escape™ Features include whisper quiet operation and dramatically reduced conducted noise levels. Also includes EPR with Easy-Breathe technology; Includes usage-based compliance downloads for accurate patient monitoring along with a standard data card for easy two-way data transfer. Optional accessories support wired and wireless compliance data transfer.</p>	E0601
36011	<p>S9 Escape & H5i S9 Escape packaged with the H5i humidifier (Item #36900).</p>	E0601, +E0562
36021	<p>S9 Escape & Climate Control Kit S9 Escape packaged with the H5i humidifier and the ClimateLine Tubing (Item #36995)</p>	E0601, +E0562, +A4604

Item #	Description	HCPCS
36004	S9 VPAP™ S A fully-integrated bilevel device, including EPR with Easy-Breathe technology, simple menu settings and color LCD. Features whisper quiet operation and dramatically reduced noise levels. AHI with Central sleep apnea (CSA) detection. Optional oximetry monitoring, wireless data compliance transfer.	E0470
36014	S9 VPAP S & H5i S9 VPAP S packaged with the H5i humidifier (Item #36900).	E0470, +E0562
36024	S9 VPAP S & Climate Control Kit S9 VPAP S with the H5i humidifier and the ClimateLine tubing (Item #36995).	E0470, +E0562 +A4604
36006	S9 VPAP™ Auto A fully-integrated auto-adjusting bilevel device with enhanced AutoSet algorithm continually monitors breathing. Simple menu settings and color LCD display. It also features a sleep quality indicator, which allows patients to view their therapy progress. Enhanced Vsync™ technology compensates for leak and provides superior patient-ventilator synchrony. TiControl™ (Ti Max/Min) feature accommodates individual respiratory conditions. AHI with central sleep apnea (CSA) detection. Optional oximetry monitoring and wireless compliance data transfer.	E0470
36016	S9 VPAP Auto & H5i S9 VPAP Auto packaged with the H5i humidifier (Item #36900).	E0470, +E0562
36026	S9 VPAP Auto & Climate Control Kit S9 VPAP Auto with the H5i humidifier and the ClimateLine Tubing (Item #36995)	E0470, +E0562 +A4604
36008	S9 VPAP™ ST A full-integrated bilevel device that delivers comfortable and effective therapy with exceptional patient-ventilator synchrony to patients requiring backup rate. It also features optional integrated oximetry and a sleep quality indicator, which allows patients to view their therapy progress on a daily basis.	E0471
36018	S9 VPAP ST & H5i S9 VPAP ST packaged with the H5i humidifier (Item #36900).	E0471 +E0562
36028	S9 VPAP ST & Climate Control Kit S9 VPAP ST with the H5i humidifier and the ClimateLine tubing (Item #36995).	E0471 +E0562 +A4604
36007	VPAP™ Adapt An adaptive servo-ventilator specifically designed to treat central sleep apnea (CSA) in all its forms, including mixed events. It also features optional integrated oximetry and a sleep quality indicator, which allows patients to view their therapy progress on a daily basis.	E0471
36017	VPAP Adapt & H5i™ VPAP Adapt packaged with the H5i humidifier (Item #36900).	E0471 +E0562
36027	VPAP Adapt & Climate Control Kit VPAP Adapt with the H5i humidifier and the ClimateLine tubing (Item#36995).	E0471 +E0562 +A4604

Item #	HUMIDIFIERS	HCPCS
36900	H5i - for use with S9 Flow Generators	E0562
36901	H5i with cleanable water tub - for use with S9 Flow Generators	E0562
H5i HUMIDIFIER ACCESSORIES		
36995	H5i ClimateLine Tubing - 1 pack	A4604
36997	H5i ClimateLine^{MAX}™ Tubing - 1 pack	A4604
36802	H5i Water Tub	A7046
36800	H5i Cleanable Water Tub	A7046
36891	H5i Flip Lid	N/A
36895	H5i for Her Flip Lid	N/A
36892	H5i Flip Lid Seal	N/A
36894	H5i Air Outlet and Outlet Clip	N/A
S9 FLOW GENERATOR ACCESSORIES		
36810	SlimLine™ Tubing	A7037
36811	SlimLine Tubing Wrap	N/A
24947	ResMed Power Station II Battery - requires S9 90W Power Supply to charge	N/A
36821	S9 90W Power Supply Unit Plus Power Cord - can be used with or without humidifier - for use with S9 and ResMed Power Station Station II Battery	N/A
36920	S9 30W Power Supply Unit	N/A
36860	S9 Travel Bag	N/A
36861	S9 for Her Travel Bag	N/A
36859	S9 Filter Cover	N/A
36970	DC/DC Converter - 24V/90W output	N/A
36955	S9 Embletta Adapter	N/A
S9 FILTERS		
36850	S9 Standard Filter - 1 pack	A7038
36851	S9 Standard Filter - 2 pack	A7038
36852	S9 Standard Filter - 12 pack	A7038
36853	S9 Standard Filter - 50 pack	A7038
36855	S9 Hypo-allergenic Filter - 1 pack	A7038
36856	S9 Hypo-allergenic Filter - 2 pack	A7038
36857	S9 Hypo-allergenic Filter - 12 pack	A7038
36858	S9 Hypo-allergenic Filter - 50 pack	A7038



Item #	Description	HCPCS
24157	Stellar™ 100 Adult and Pediatric Ventilator Ultra-premium non-invasive ventilator with invasive option for adult and pediatric (13 kg) ventilation. Pressure from 3–40 cm H ₂ O, PS to 37 cm H ₂ O, and internal battery. FiO ₂ and SpO ₂ monitoring (optional)	E0471

Item #	Stellar Humidifier and Accessories	HCPCS
24927	H4i Heated Humidifier for Stellar	E0562
36810	SlimLine tube - for use with S9 and Stellar	A7037
24977	Stellar Mobility bag	N/A
22305	XPOD Oximeter	N/A
16936	Oximeter Sensor Cable	N/A
24955	FiO₂ Monitoring Kit	N/A
CAP008617	FiO₂ Monitoring Sensor	N/A
24976	ResMed Leak Port	N/A
24934	Stellar Filter - hypoallergen - 12 pack	A7038
24905	Standard stand (without top plate)	N/A
24913	VPAP Series top plate	N/A
24917	Premium circuit support arm	N/A
24919	Tube clip–double	N/A
24920	Tube clip–single	N/A

Item #	Description	HCPCS
24116	S7 VPAP™ III ST-A Premium non-invasive ventilator with expanded pressure range from 3–30 cm H ₂ O and a comprehensive suite of alarms.	E0471
36008	S9 VPAP™ ST A full-integrated bilevel device that delivers comfortable and effective therapy with exceptional patient-ventilator synchrony to patients requiring backup rate. It also features optional integrated oximetry and a sleep quality indicator, which allows patients to view their therapy progress on a daily basis.	E0471
36018	S9 VPAP ST & H5i S9 VPAP ST packaged with the H5i humidifier (Item #36900).	E0471 +E0562
36028	S9 VPAP ST & Climate Control Kit S9 VPAP ST with the H5i humidifier and the ClimateLine tubing (Item #36995).	E0471 +E0562 +A4604

Item #	C-SERIES FLOW GENERATOR	HCPCS
34000	C-Series Tango™ Features include classic CPAP functionality with ResMed superior quality and reliability. Optional humidification and disposable water chamber.	E0601



Item #	COMPLIANCE SOFTWARE	HCPCS
31313	ResScan software - clinical management software for use with all ResMed flow generator products.	N/A
S9 DATA ACCESSORIES		
36840	Protective Folder with SD card - 10 pack	A9279
36845	Protective Folder without SD card - 10 pack	N/A
36931	SD Card Reader	N/A
S8 DATA ACCESSORIES		
22903	ResScan data cards with protective folder - 10 pack	A9279
22906	Protective Folder (no data card) - 10 pack	N/A
22207	ResScan card reader	N/A
22212	SmartMedia card reader II - USB	N/A
30942	SmartMedia card - 2 pack	A9279
30947	SmartMedia card spare labels - 10 pack	N/A
22201	Serial port adapter	N/A
22203	ResScan USB serial adapter	N/A
22915	Data module cover	N/A
17952	Null modem cable - for use with VPAP III, AutoSet Spirit™, AutoSet Respond™, S7™ Elite CPAP, S8 CPAP and AutoSet T flow generators	N/A
S9 MODULE AND OXIMETRY		
36930	S9 Wireless Module	A9279
36950	S9 USB Adapter	N/A
36940	S9 Oximeter Adapter	N/A
369100	S9 Complete Oximetry Kit (includes Item #36940 and Item #36942)	N/A
S8 MODULE AND OXIMETRY		
22206	S8 ResLink™ Includes 32 MB SmartMedia card, user's manual and mailing bag.	A9279
22230	S8 ResTraxx GSM device - Includes a carry bag and Quick Reference Guide Card. A proprietary system comprised of both hardware and software that facilitates the wireless transfer of patient compliance data via cellular technology to an internet-accessed portal known as ResTraxx™ Online. There is also a monthly subscription fee.	A9279
222102	S8 Complete Oximetry Kit (includes Item #22206 and Item #1430050)	A9279
S7 MODULE AND OXIMETRY		
30924	S7 ResLink PMM - patient management module Includes 32 MB SmartMedia card, user's manual and mailing bag.	A9279
309100	S7 Complete Oximetry Kit (includes Item #30924 and Item #1430050)	A9279



Item #	OXIMETRY ACCESSORIES	HCPCS
1430050	Oximeter Spirit PMM kit ResLink cable and XPod for S7& S8	N/A
36942	S9 Oximetry Accessory Kit (includes reusable soft sensor, 3M cable, sensor clip, Xpod)	N/A
707560	Reusable soft sensor - 3M	N/A
1431002	Disposable flex sensor - 3M - includes sensor and item #70276 (adhesive strips, single use - 25 pk)	N/A
22961	ResLink converter cable - This converter cable is for use with ApneaLink™ XPOD to connect the oximeter port for either the S7 ResLink or S8 ResLink modules.	N/A
36941	S9 Oximeter converter cable - This converter cable is for use with existing ResLink XPOD to connect with the S9 Oximeter module.	N/A



Item #	Description	HCPCS
22328	<p>ApneaLink™ Plus Complete Kit A portable home sleep test diagnostic device that records up to four channels of information: respiratory effort, pulse, oxygen saturation and nasal flow. Kit includes ApneaLink Plus recorder, software CD, quick software setup guide, carry bag, EasySense respiratory effort sensor, 2 reusable belts, USB download cable, 2 AA batteries, 3 nasal cannulas, XPOD oximeter, 3 single use oximetry sensors and sensor clip.</p>	N/A
22319	<p>ApneaLink Plus Basic Kit A portable home sleep test diagnostic device that records up to four channels of information: respiratory effort, pulse, oxygen saturation and nasal flow. Kit includes ApneaLink Plus recorder, software CD, quick software setup guide, carry bag, EasySense respiratory effort sensor, 2 reusable belts, USB download cable, 2 AA batteries and 1 nasal cannula.</p>	N/A
22302	<p>ApneaLink™ A sleep screening device with up to three channels (with oximetry) used to objectively assist in identifying the risk for sleep-disordered breathing. Includes ApneaLink software CD, operating instructions, carry bag, 3 nasal cannulas, single use belt, USB download cable and 2 AA batteries.</p>	N/A
22304	<p>ApneaLink Oximetry Kit Provides oximetry and pulse information to the ApneaLink/ApneaLink Plus recorder devices. Uses single-use or reusable probes. Must be used in conjunction with ApneaLink and ApneaLink Plus only. Kit includes oximeter (XPOD), oximeter sensor clip and 3 single use sensors.</p>	N/A

Item #	ACCESSORIES AND CONSUMMABLES	HCPCS
22336	ApneaLink Plus and ApneaLink Software CD 9.2 (compatible with all devices)	N/A
22333	ApneaLink Plus EasySense Respiratory Effort Sensor (for use with ApneaLink Plus ONLY)	N/A
NASAL CANNULAS		
70388	ApneaLink standard cannula - 25 pack	N/A
70319	ApneaLink nasal/oxygen cannula - 25 pack	N/A
BELTS		
629052	Reusable belt - for use with ApneaLink and ApneaLink Plus	N/A
70406	Single use belt - 10 pack (for use with ApneaLink ONLY)	N/A
MISCELLANEOUS		
629051	USB download cable	N/A
22309	Rechargeable battery kit - includes battery charger and 4 AA rechargeable batteries	N/A
OXIMETRY COMPONENTS		
22306	Oximeter sensor clip	N/A
70413	Reusable soft sensor - large	N/A
70567	Reusable soft sensor - small	N/A
70726	Reusable flex sensor kit, 36" includes sensor and Item #70276	N/A
70276	Adhesive strips, single use - for use with Item #1431002 or Item #70726 - 25 pack	N/A
22338	Single use sensor, adult, 36" (1 meter) (ResMed Standard) - 24 pack. Replaces Item #22316	N/A
22337	Single use sensor, 36" (1 meter) (value line) - 24 pack. Replaces Item #70412	N/A

Flow Generator Accessories and Replacement Parts

Item #	POWER SUPPLIES	HCPCS
70253	110V inverter for Sullivan™ CPAP and VPAP systems	N/A
70199	Battery adapter, 12 volt	N/A
26932	DC converter - for VPAP ST, VPAP Auto 25 and VPAP S	N/A
22015	DC converter - for VPAP III ST-A	N/A
33942	DC-12 converter - for S8 CPAP unit and VPAP Auto	N/A
33938	Power cord - for S8 CPAP and C-Series Tango units	N/A
14951	Power cord - for Sullivan V and ResMed S6 CPAP units, S7 Elite CPAP, VPAP, AutoSet T, AutoSet Spirit and AutoSet Respond units	N/A
26905	Power cord - for VPAP Adapt SV and VPAP III ST-A	N/A
33948	S8 connector plug - 10 pack	N/A
Item #	RESMED POWER STATION II	HCPCS
24947	ResMed Power Station II Battery - requires S9 90W Power Supply to charge	N/A
24923	ResMed Power Station II - Battery and bag	N/A
24959	DC cable, Stellar	N/A
24981	DC cable, Stellar, 0.5m	N/A
24961	DC cable, S9 Series	N/A
24960	DC cable, S8 VPAP	N/A
24946	PSU adaptor, ResMed Power Station II	N/A
36821	90W Power Supply Unit Plus Power Cord - can be used with or without humidifier - for use with S9 and ResMed Power Station II Battery	N/A
24964	Battery coupler kit, ResMed Power Station II	N/A
24965	Battery coupler, ResMed Power Station II	N/A
24952	Carry bag, ResMed Power Station II	N/A
24943	Velcro strap, ResMed Power Station II	N/A
TUBINGS AND AIR DELIVERY SYSTEMS		
26923	Air delivery system - for VPAP Adapt SV Includes sensor tubing, tube clips, proximal cuff and standard 2 meter air hose	N/A
14911	Air tubing connector - 5 pack	N/A
14994	Air tubing, standard - light gray, ribbed (6')	A7037
14999	Air tubing, clear gray ribbed (9', 2.75 meter)- for use with all S7, S8 products and S9 VPAP only	A7037
14986	Air tubing, autoclave - dark gray (6'6", 2 meter)	A7037
14922	Air tubing, gray (9'9", 3 meter)	A7037
14987	Air tubing, standard - light gray, cuffed (6'6", 2 meter)	A7037
26960	Oxygen adapter port - for VPAP ST, VPAP Auto 25 and VPAP S	N/A
14915	Oxygen sideport connector	N/A
26931	Proximal cuff - for VPAP Adapt SV	N/A
26924	Sensor tubing with Luer lock - for VPAP Adapt SV	N/A
26925	Tube clips - 9 pack - for VPAP Adapt SV	N/A
CARRY / TRAVEL BAGS		
33925	Carry pouch for S8	N/A
34902	C-Series Tango carry bag	N/A
33972	S8 Series II premium travel bag (black)	N/A
33956	S8 standard travel bag (blue)	N/A
30923	VPAP Adapt SV carry bag	N/A

AIR FILTER COVERS

21937	ResMed S6 air filter cover	N/A
30909	S7 Series and VPAP III air filter cover	N/A
33923	Filter cover for S8	N/A
33959	Blue filter cover for S8	N/A

MANOMETERS AND PRESSURE SENSOR ADAPTERS

14535	Manometer, gauge type and portable kit Includes carrying case, tubing and stand. Range: 0–25 cm H ₂ O, 0.5 increments	N/A
30901	Pressure sensor adapter - clear plastic adapter for external humidifier connection with the AutoSet Spirit, AutoSet Respond and S7 Elite	N/A
17947	Pressure sensor adapter - for use with AutoSet T	N/A
30916	Pressure sensor adapter - rubber adapter for AutoSet Spirit, AutoSet Respond and S7 Elite CPAP	N/A

Item #	Description	HCPCS
33915	S8, S8 II and C-Series Tango filters - 1 pack	A7038
33973	S8, S8 II and C-Series Tango filters - 2 pack	A7038
33974	S8, S8 II and C-Series Tango filters - 12 pack	A7038
33918	S8, S8 II and C-Series Tango filters - 50 pack	A7038
33919	S8, S8 II and C-Series Tango hypo filters - 1 pack	A7038
33975	S8, S8 II and C-Series Tango hypo filters - 2 pack	A7038
33976	S8, S8 II and C-Series Tango hypo filters - 12 pack	A7038
33922	S8, S8 II and C-Series Tango hypo filters - 50 pack	A7038
30919	S7 Series and VPAP III filters - 1 pack	A7038
30907	S7 Series and VPAP III filters - 3 pack	A7038
30908	S7 Series and VPAP III filters - 50 pack	A7038
21944	ResMed S6 air filter strips - 1 pack	A7038
21935	ResMed S6 air filter strips - 3 pack	A7038
21936	ResMed S6 air filter strips - 50 pack	A7038
21946	Air filter strips - fine, hypoallergenic - 1 pack	A7038
21926	Air filter strips - fine, hypoallergenic, fits all Sullivan V CPAP and VPAP III units - 5 pack	A7038
21947	Air filter strips - fine, hypoallergenic, fits all ResMed S6 units - 10 pack	A7038
17957	AutoSet T filter - 1 pack	A7038
17937	AutoSet T filter - 3 pack	A7038
17938	AutoSet T filter - 50 pack	A7038
14985	Air filter strips fits all Sullivan V, CPAP and VPAP II units - 1 pack	A7038
14907	Air filter strips fits all Sullivan V, CPAP and VPAP II units - 50 pack	A7038
14908	Air filter strips fits all Sullivan V, CPAP and VPAP II units - 3 pack	A7038
14909	Air filter assembly - fits all Sullivan V, CPAP and VPAP II units	N/A
14910	Air filter assembly - fine, hypoallergenic, fits all Sullivan V, CPAP and VPAP II units - 1 pack	N/A
26926	VPAP Adapt SV air inlet filter, standard - 3 pack	N/A
26927	VPAP Adapt SV air inlet filter, hypoallergenic - 3 pack	N/A





Item #	Description	HCPCS
26940	H4i - for use with all S8, S8 II and VPAP bilevel products.	E0562
26950	H4i with optional cleanable water chamber - for use with all S8, S8 II and VPAP bilevel products.	E0562
30902	HumidAire 2i™ - for use with S7™ Lightweight CPAP, S7 Elite CPAP, AutoSet Spirit, AutoSet Respond, VPAP III and VPAP Adapt SV	E0562
34900	C-Series heated humidifier Features disposable water chamber	E0562
Item #	Description	HCPCS
26952	H4i water chamber (standard)	A7046
26958	H4i water chamber (cleanable)	A7046
26959	H4i water chamber upgrade kit	A7046
26953	H4i top cover and seal	N/A
26954	H4i tub plate seal	N/A
26955	H4i cleaning plug	N/A
26963	H4i spill guard	N/A
33926	HumidAire 3i top cover and seal	N/A
33927	HumidAire 3i top cover seal	N/A
33935	HumidAire 3i carry pouch	N/A
33947	HumidAire 3i cleaning plug pack	N/A
30951	HumidAire 2i tub assembly	A7046
30915	HumidAire 2i water chamber o-ring	N/A
30914	HumidAire 2i water chamber tub plate	N/A
30938	HumidAire 2i water chamber cross member, PC	N/A
30961	HumidAire 2i base PC-side catch	N/A
30939	HumidAire 2i latch, PC	N/A
30955	HumidAire 2i side catches, PC	N/A
30956	HumidAire 2iC tub assembly	A7046
30912	HumidAire 2i and 2iC water chamber seal	N/A
34901	C-Series water chamber	A7046
34903	C-Series adapter seal - 10 pack	N/A
18920	HumidAire water chamber	A7046
18913	HumidAire water chamber base	N/A
18911	HumidAire water chamber lid	N/A
18912	HumidAire water chamber seal	N/A
18910	Passover humidifier connector	N/A
18919	Short air tubing - 19 cm/7.5 inch and white adapter	N/A
33963	Tubing wrap (6' 6") - reduces rainout in tubing caused by heated humidification	N/A
18916	Tubing (20 inch) - use with HumidAire humidifier	A7037

This Home Care Provider Product List has been prepared solely for the use of ResMed's customers.

ResMed reserves the right to modify any of the information contained herein at any time without notice, including pricing, product descriptions and configurations, and warranty terms.

All information contained in this Product List is confidential. You agree not to disclose such information to any third party without the prior written approval of ResMed. While ResMed has taken reasonable precautions to ensure the accuracy of the pricing and information presented in this document as of the date of publication, ResMed strongly recommends that you contact customer service and/or a sales representative to confirm the accuracy of the information or to see if other discounts, promotions, or specials may apply.

All purchases of ResMed products by Customer will be governed by ResMed's Standard Terms and Conditions of Sales located on the back of each invoice, unless otherwise agreed upon in separate writing by ResMed.

Please contact ResMed's Customer Service Department at (800) 424-0737 if you have any questions or need information about ResMed's standard terms and conditions.



ResMed Corp San Diego, CA, USA +1 858 836 5000 or 1 800 424 0737 (toll free), ResMed Ltd Bella Vista, NSW, Australia +61 (2) 8884 1000 or 1 800 658 189 (toll free). Offices in Austria, Brazil, Finland, France, Germany, Hong Kong, Japan, Malaysia, Mexico, Netherlands, New Zealand, Norway, Singapore, Spain, Sweden, Switzerland, United Kingdom (see website for details). Protected

Activa, AutoScan, AutoSet, AutoSet Respond, AutoSet Spirit, ClimateLine, ConvertAble, EPR, Gecko, H4i, HumidAire, Kidsta, Liberty, Mirage, Mirage Activa, Mirage FX, Mirage Kidsta, Mirage Liberty, Mirage Micro, Mirage Quattro, Mirage Swift, Mirage Vista, Pixi, Quattro, Rescontrol, ResLink, ResScan, ResTraxx, SmartStart, S7 Elite, S7 Lightweight, S8, S8 Compact, S8 Elite, S8 Escape, S9, S9 AutoSet, S9 Elite, S9 Escape, S9 Escape Auto, SlimLine, Stellar, Sullivan, Swift, Tango, Ultra Mirage, Vista, VPAP, VPAP Adapt SV, VPAP Auto, VPAP Malibu and VSet are trademarks of ResMed Ltd.

Activa, ApneaLink, AutoSet, Escape, HumidAire, Kidsta, Mirage, Mirage Quattro, SmartStart, Sullivan, Vista and VPAP are registered in U.S. Patent and Trademark Office.

ApneaLink is a trademark of MAP Medizin Technologie GmbH and Breath-O-Prene™ is a trademark of AccuMed. ©2007 ResMed Ltd.